



ADSS Antibody

Product Code	CSB-PA001398GA01HU
Abbreviation	ADSS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30520
Immunogen	Human ADSS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	adenylosuccinate synthase;ADSS;ADEH;MGC20404 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ADSS
Target Details	Adenylosuccinate synthetase catalyzes the first committed step in the conversion of IMP to AMP