



Recombinant *Oryza sativa* subsp. *indica* Actin-related protein 2 (ARP2)

Product Code	CSB-EP001248OFF-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2YUL5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDSGNVVVCD NGTGYVKCGF AGENFPTSVF PCVVGRPLLR YEESLQEDEL TDIVVGAACA DLRHQLDVSY PVTNGIVQSW DDMGHIWDHA FYSELKVDPS ECKILLTDPD LNPVKNREKM IETMFEKYNF AGVFIQVQAV LSLYAQGLLT GLVIDSGDGV THVVPVVDGF SYPHITKRMN VAGRHITSYL VDLLSRRGYA MNKSADFETV REIKEKLCYL SYDYKREYQL GLETTILVKS YTLDPGRVIK VGTERFQAPE ALFTPELIDV EGDGMADMAF RCIQEMDIDN RMTLYQHIVL SGGSTMY PGL PSRLEKEMLD RYLDVVLKGN KDGLKCLRRLR IEDPPRRKHM VYLGGA VLAG IMKDAPEFWI TRQEYQEEGL ACLRKCGQA
Source	E.coli
Target Names	ARP2
Protein Names	Recommended name: Actin-related protein 2
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.