



AGFG1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA001442GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P52594 |
| Immunogen | Human AGFG1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | ArfGAP with FG repeats 1;AGFG1;DKFZp686I15205;HRB;MGC116938;MGC116940;RAB;RIP ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | AGFG1 |
| Target Details | This protein is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene. |