



ALG1 Antibody

Product Code	CSB-PA001593GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT22
Immunogen	Human ALG1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	asparagine-linked glycosylation 1, beta-1,4-mannosyltransferase homolog (S. cerevisiae);ALG1;HMAT1;HMT-1;HMT1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ALG1
Target Details	The enzyme encoded by this gene catalyzes the first mannosylation step in the biosynthesis of lipid-linked oligosaccharides. This gene is mutated in congenital disorder of glycosylation type I _k .
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.