



ALG5 Antibody

Product Code	CSB-PA001602GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y673
Immunogen	Human ALG5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	asparagine-linked glycosylation 5, dolichyl-phosphate beta-glucosyltransferase homolog (<i>S. cerevisiae</i>);ALG5;RP11-421P11.2;bA421P11.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ALG5
Target Details	This gene encodes a member of the glycosyltransferase 2 family. The encoded protein participates in glycosylation of the oligomannose core in N-linked glycosylation of proteins. The addition of glucose residues to the oligomannose core is necessary to ensure substrate recognition, and therefore, effectual transfer of the oligomannose core to the nascent glycoproteins. Multiple transcript variants encoding different isoforms have been found for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.