



Recombinant Danio rerio Class A basic helix-loop-helix protein 9 (bhlha9)

Product Code	CSB-MP002685DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3KNX5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MTPVSVCRFS MASRGSFTGS EFSEEEP DGS LQSELDSSD GLGGLNEPEE RLIKKRSRPV RSKARRVAAN VRERKRILDY NQAFNALRVA LHHDLGKRL SKIATLQRAI NRISALSVFL TNNPPVGVAK PCGHLECCQPG GLWAEAEPSM QNFLTWHQPL NQHLQTSIHR LSSEQHVFTG PACPPSPHYP CFSPDTQLYP AASVSPPPRY GRIGDVGAYQ QGVWGGNHAD GYGESLQTLF LPWHMGYLQD AGLQHCPNTL CYTQCKTVRA L
Source	Mammalian cell
Target Names	bhlha9
Protein Names	Recommended name: Class A basic helix-loop-helix protein 9 Short name= bHLHa9 Alternative name(s): Class B basic helix-loop-helix factor 42 Short name= bHLHf42
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.