



ANKRD2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA001751GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9GZV1 |
| Immunogen | Human ANKRD2 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Form | Liquid |
| Storage Buffer | PBS with 0.02% sodium azide and 50% glycerol, pH 7.3. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ANKRD2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |