




ANXA11 Antibody

Product Code	CSB-PA001838GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50995
Immunogen	Recombinant Human Annexin A11 protein (1-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA11

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate

All lanes: ANXA11 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 54 kDa

Observed band size: 54 kDa

Exposure time: 5s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.