



AOX1 Antibody

Product Code	CSB-PA001858GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06278
Immunogen	Human AOX1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	aldehyde oxidase 1;AOX1;AO;AOH1 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AOX1
Target Details	Aldehyde oxidase produces hydrogen peroxide and, under certain conditions, can catalyze the formation of superoxide. Aldehyde oxidase is a candidate gene for amyotrophic lateral sclerosis.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.