



AP2B1 Antibody

Product Code	CSB-PA001871GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63010
Immunogen	Human AP2B1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	adaptor-related protein complex 2, beta 1 subunit;AP2B1;ADTB2;AP105B;AP2-BETA;CLAPB1;DKFZp781K0743 ;
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AP2B1
Target Details	This protein is one of two large chain components of the assembly protein complex 2, which serves to link clathrin to receptors in coated vesicles. The encoded protein is found on the cytoplasmic face of coated vesicles in the plasma membrane. Two transcript variants encoding different isoforms have been found for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.