



APBB2 Antibody

Product Code	CSB-PA001892GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92870
Immunogen	Human APBB2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	amyloid beta (A4) precursor protein-binding, family B, member 2;APBB2;DKFZp434E033;FE65L;FE65L1;MGC35575 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APBB2
Target Details	This protein interacts with the cytoplasmic domains of amyloid beta (A4) precursor protein and amyloid beta (A4) precursor-like protein 2. This protein contains two phosphotyrosine binding (PTB) domains, which are thought to function in signal transduction. Polymorphisms in this gene have been associated with Alzheimer s disease. Alternative splicing results in multiple transcript variants.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.