



APIP Antibody

Product Code	CSB-PA001906GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GX9
Immunogen	Human APIP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	APAF1 interacting protein;APIP;APIP2;CGI-29;CGI29;MMRP19;dJ179L10.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.