



Recombinant Prochlorococcus marinus Queuine tRNA-ribosyltransferase (tgt)

Product Code	CSB-MP019146PZE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3PAZ2
Product Type	Recombinant Protein
Immunogen Species	Prochlorococcus marinus (strain MIT 9301)
Purity	>85% (SDS-PAGE)
Sequence	MFEFEITSDC INTGARTGIF HTPNGKVNTP KFMPVGTLAT VKGISSKQLI STGSEMILSN TFHLHLQPGE KLVKASGGIH KFMNWPKPIL TDSGGYQVFS LAKLNNISDE GVEFKNPRDG SHVFLSPEKV IKIQMDLGS D VAMAFDHCPP HTANENDIED SLERTHSWLQ KCVETHQKSN QALFGIVQGG KYPRLREYSA KYTSSFDLPG IAVGGVSVGE AVEEHSVIV YVPKFLPINK PRYLMGIGSL KEISLAVANG FDIFDCVLPT RLGRHGTAFF NDERLNLRNA RFKNDFSPID KTCKCETCKS YSRAYLHHLI RNDEILGLSL ISLHNIAHLI RFTNAISTAI RDNCFTNDFA PWKTSSIAHH TW
Source	Mammalian cell
Target Names	tgt
Protein Names	Recommended name: Queuine tRNA-ribosyltransferase EC= 2.4.2.29 Alternative name(s): Guanine insertion enzyme tRNA-guanine transglycosylase
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.