



ARHGEF7 Antibody

Product Code	CSB-PA002056GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14155
Immunogen	Human ARHGEF7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Rho guanine nucleotide exchange factor (GEF) 7;ARHGEF7;BETA-PIX;COOL1;DKFZp686C12170;DKFZp761K1021;KIAA0142;KIAA0412;Nbla10314;P50;P50BP;P85;P85COOL1;P85SPR;PAK3;PIXB ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGEF7
Target Details	Rho GTPases play a fundamental role in numerous cellular processes triggered by extracellular stimuli that work through G protein coupled receptors. The encoded protein belongs to a family of cytoplasmic proteins that activate the Ras-like family of Rho proteins by exchanging bound GDP for GTP. It forms a complex with the small GTP binding protein Rac1 and recruits Rac1 to membrane ruffles and to focal adhesions. This protein can induce membrane ruffling. Multiple alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.