



ARL13B Antibody

Product Code	CSB-PA002074GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SXY8
Immunogen	Human ARL13B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ADP-ribosylation factor-like 13B;ARL13B;ARL2L1;DKFZp686E2075;DKFZp686L2472;DKFZp686M2074;DKFZp761H079;JBTS8;MGC120611;MGC120612 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARL13B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.