



ASCC1 Antibody

Product Code	CSB-PA002196GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9N2
Immunogen	Human ASCC1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	activating signal cointegrator 1 complex subunit 1;ASCC1;ASC1p50;CGI-18;RP11-150D20.4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ASCC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.