



ASCC3 Antibody

Product Code	CSB-PA002198GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3C0
Immunogen	Human ASCC3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	activating signal cointegrator 1 complex subunit 3;ASCC3;ASC1p200;DJ467N11.1;HELIC1;MGC26074;RNAH;dJ121G13.4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ASCC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.