



ATG3 Antibody

Product Code	CSB-PA002288GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NT62
Immunogen	Human ATG3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ATG3 autophagy related 3 homolog (S. cerevisiae);ATG3;APG3;APG3-LIKE;APG3L;DKFZp564M1178;FLJ22125;MGC15201;PC3-96 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.