



ATOH1 Antibody

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA002306GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q92858 |
| Immunogen | Human ATOH1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | atonal homolog 1 (Drosophila);ATOH1;ATH1;HATH1;MATH-1;bHLHa14 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ATOH1 |
| Target Details | This protein belongs to the basic helix-loop-helix (BHLH) family of transcription factors. It activates E-box dependent transcription along with E47. |