



Recombinant Mannose-1-phosphate guanyltransferase beta (tag-335)

Product Code	CSB-BP009579CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3QMC8
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MKALILVGGY GTRLRPLTLT QPKPLVEFAN KPMMLHQMEA LAEVGVDTVV LAVSYRAEQL EQEMTVHADR LGVKLIFSLE EEPLGTAGPL ALARKHLEGD APFFVLNSDV ICDFFPKQMV EFHKNHGKEG TIAVTKVEEP SKYGVVVFQD DKGKIDDFVE KPQEYVGNKI NAGLYIFSSK ILDRIPKPT SIEKEIFPEM AFSGNLYAFV LPGFWMDVGQ PKDFLKGMSL FLNHCHTTKS DKLETGSNIH PTATIRGNVM VDPSATVGEN CVIGPDVVIG PRVKIEGGVR ILHSTILSDS SIGNYSWVSG SIVGRKCHIG SWVRIENICV IGDDVVVKDE LYLNGASVLP HKSIAVNVPS KDIIM
Source	Baculovirus
Target Names	tag-335
Protein Names	Recommended name: Mannose-1-phosphate guanyltransferase beta EC= 2.7.7.13 Alternative name(s): GDP-mannose pyrophosphorylase B GTP- mannose-1-phosphate guanyltransferase beta
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.