



# Recombinant Mannose-1-phosphate guanyltransferase beta (tag-335)

<b>Product Code</b>	CSB-EP009579CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A3QMC8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKALILVGGY GTRLRPLTLT QPKPLVEFAN KPMMLHQMEA LAEVGVDTVV LAVSYRAEQL EQEMTVHADR LGVKLIFSLE EEPLGTAGPL ALARKHLEGD APFFVLNSDV ICDFPFKQMV EFHKNHGKEG TIAVTKVEEP SKYGVVVFQD DKGKIDDFVE KPQEYVGNKI NAGLYIFSSK ILDRIPKPT SIEKEIFPEM AFSGNLYAFV LPGFWMDVGQ PKDFLKGMSL FLNHCHTTKS DKLETGSNIH PTATIRGNVM VDPSATVGEN CVIGPDVVIG PRVKIEGGVR ILHSTILSDS SIGNYSWVSG SIVGRKCHIG SWVRIENICV IGDDVVVKDE LYLNGASVLP HKSIAVNVPS KDIIM
<b>Source</b>	E.coli
<b>Target Names</b>	tag-335
<b>Protein Names</b>	Recommended name: Mannose-1-phosphate guanyltransferase beta EC= 2.7.7.13 Alternative name(s): GDP-mannose pyrophosphorylase B GTP- mannose-1-phosphate guanyltransferase beta
<b>Expression Region</b>	1-365
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.