



Recombinant Mouse Apolipoprotein D (Apod)

Product Code	CSB-MP001935MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51910
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	QNFHLGKCPS PPVQENFDVK KYLGRWYEIE KIPASFEKGN CIQANYSLME NGNIEVLNKE LSPDGTMNQV KGEAKQSNVS EPAKLEVQFF PLMPPAPYWI LATDYENYAL VYSCTTFFWL FHVDFVWILG RNPYLPPETI TYLKDILTSN GIDIEKMTT DQANCPDFL
Source	Mammalian cell
Target Names	Apod
Protein Names	Recommended name: Apolipoprotein D Short name= Apo-D Short name= ApoD
Expression Region	21-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a component of high density lipoprotein that has no marked similarity to other apolipoprotein sequences. It has a high degree of homology to plasma retinol-binding protein and other members of the alpha 2 microglobulin protein superfamily of carrier proteins, also known as lipocalins. This glycoprotein is closely associated with the enzyme lecithin:cholesterol acyltransferase - an enzyme involved in lipoprotein metabolism.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.