



B3GALNT2 Antibody

Product Code	CSB-PA002489GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCR0
Immunogen	Human B3GALNT2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	beta-1,3-N-acetylgalactosaminyltransferase 2;B3GALNT2;B3GalNAc-T2;MGC39558 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	B3GALNT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.