



Recombinant Human Ubiquitin-conjugating enzyme E2 E1 (UBE2E1)

Product Code	CSB-BP025448HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P51965
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SDDDSRAST SSSSSSSSNQ QTEKETNTPK KKEKVSMSK NSKLLSTSAK RIQKELADIT LDPPPNCASAG PKGDNIYEWR STILGPPGSV YEGGVFFLDI TFTPEYPFKP PKVTFRTRIV HCNINSQGVI CLDILKDNWS PALTISKVLL SICSLLTDCN PADPLVGSIA TQYMTNRAEH DRMARQWTKR YAT
Source	Baculovirus
Target Names	UBE2E1
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 E1 EC= 6.3.2.19 Alternative name(s): UbchH6 Ubiquitin carrier protein E1 Ubiquitin-protein ligase E1
Expression Region	2-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.