



BCAT2 Antibody

Product Code	CSB-PA002601GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15382
Immunogen	Human BCAT2
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	branched chain aminotransferase 2, mitochondrial;BCAT2;BCAM;BCT2 ;
Immunogen Species	Homo sapiens (Human)
Target Names	BCAT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.