



BHLHE40 Antibody

Product Code	CSB-PA002691GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14503
Immunogen	Human BHLHE40
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	basic helix-loop-helix family, member e40;BHLHE40;BHLHB2;DEC1;FLJ99214;SHARP-2;STRA13;Stra14 ;SHARP2
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BHLHE40
Target Details	This gene encodes a basic helix-loop-helix protein expressed in various tissues. Expression in the chondrocytes is responsive to the addition of Bt2cAMP. The encoded protein is believed to be involved in the control of cell differentiation.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.