



BLVRB Antibody

Product Code	CSB-PA002722GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30043
Immunogen	Human BLVRB
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	biliverdin reductase B (flavin reductase (NADPH));BLVRB;BVRB;FLR;MGC117413;SDR43U1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BLVRB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.