



Recombinant Human WD repeat-containing protein 91 (WDR91)

Product Code	CSB-YP026077HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4D1P6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MAEAVERTDE LVREYLLFRG FTHTLRQLDA EIKADKEKGF RVDKIVDQLQ QLMQVYDLAA LRDYWSYLER RLFSRLEDIY RPTIHKLKTS LFRFYLVYTI QTNRNDKAQE FFAKQATELQ NQAEWKDWFV LPFLPSPDTN PTFATYFSRQ WADTFIVSLH NFLSVLFQCM PVPVILNFDA ECQRNQQVQE ENEVLRQKLF ALQAEIHRK KEEQQPEEEE ALVQHKLPY VSNMDRLGDS ELAMVCSQRN ASLSQSPRVG FLSSLLPQSK KSPSRLSPAQ GPPQPQSSAK KESFGGQGTK GKDPTSGAKD GKSLLSGLAT GESGWSQHRQ RRLQDHGKER KELFSTTTTSQ CAEKKPEASG PEAEPCEPLH TEPVEPLTRA SSAGPEGGGV RPEQPFIVLG QEEYGEHHSS IMHCRVDCSG RRVASLDVDG VIKVVSFNPI MQTKASSISK SPLLSLEWAT KRDRLLLLGS GVGTVRLYDT EAKKNLCEIN INDNMPRILS LACSPNGASF VCSAAAPSLT SQVDFSAPDI GSKGMNQVPG RLLLWDTKTM KQQLQFSLDP EPIAINCTAF NHNGNLLVTG AADGVIRLFD MQQHECAMSW RAHYGEVYSV EFSYDENTVY SIGEDGKFIQ WNIHKSLKV SEYSLPSDAT GPFVLSGYSG YKQVQVPRGR LFAFDSEGNV MLTCSATGGV IYKLGDEKV LESCLSLGGH RAPVVTVDWS TAMDCGTCLT ASMDGKIKLT TLLAHKA</p>
Source	Yeast
Target Names	WDR91
Protein Names	Recommended name: WD repeat-containing protein 91
Expression Region	1-747
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.