



BTN3A3 Antibody

Product Code	CSB-PA002875GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00478
Immunogen	Human BTN3A3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	butyrophilin, subfamily 3, member A3;BTN3A3;BTF3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BTN3A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.