



ANAPC15 Antibody

Product Code	CSB-PA003014GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60006
Immunogen	Human C11orf51
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 11 open reading frame 51;C11orf51;DKFZp564M082;HSPC020 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ANAPC15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.