



C14orf166 Antibody

Product Code	CSB-PA003165GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y224
Immunogen	Human C14orf166
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 14 open reading frame 166;C14orf166;CGI-99;CGI99;CLE;CLE7;LCRP369;RLLM1 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RTRAF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.