



FAM195A Antibody

Product Code	CSB-PA003276GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUT9
Immunogen	Human C16orf14
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 16 open reading frame 14;C16orf14;MGC104301;MGC15416;c349E10.1 ;; FAM195A
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MCRIP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.