



SARG Antibody

Product Code	CSB-PA003507GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW04
Immunogen	Human C1orf116
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 1 open reading frame 116;C1orf116;DKFZp666H2010;FLJ36507;MGC2742;MGC4309;SARG ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SARG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.