



C1QTNF3 Antibody

Product Code	CSB-PA003648GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXJ4
Immunogen	Human C1QTNF3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	C1q and tumor necrosis factor related protein 3;C1QTNF3;C1ATNF3;CORS26;CTRP3;Corcs;Cors;FLJ37576 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	C1QTNF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.