



RTCB Antibody

Product Code	CSB-PA003810GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3I0
Immunogen	Human C22orf28
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromosome 22 open reading frame 28;C22orf28;DJ149A16.6;HSPC117;RP1-149A16.6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RTCB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.