



CACNA1B Antibody

Product Code	CSB-PA004398GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00975
Immunogen	Human CACNA1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	calcium channel, voltage-dependent, N type, alpha 1B subunit;CACNA1B;BIII;CACNL1A5;CACNN;Cav2.2 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CACNA1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.