



CALM1 Antibody

Product Code	CSB-PA004445GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62158
Immunogen	Human CALM1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	calmodulin 1 (phosphorylase kinase, delta);CALM1;CALML2;CAMI;DD132;PHKD ;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CALM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.