



# Recombinant Human DNA-directed RNA polymerases I, II, and III subunit RPABC3 (POLR2H)

<b>Product Code</b>	CSB-BP018335HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P52434
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AGILFEDIF DVKDIDPEGK KFDRVSR LHC ESEFKMDLI LDVNIQIYPV DLGDKFRLVI ASTLYEDGTL DDGEYNPTDD RPSRADQFEY VMYGKVYRIE GDETSTEAAAT RLSAYVSYGG LLMRLQGDAN NLHGFEVDSR VYLLMKKLAF
<b>Source</b>	Baculovirus
<b>Target Names</b>	POLR2H
<b>Protein Names</b>	Recommended name: DNA-directed RNA polymerases I, II, and III subunit RPABC3 Short name= RNA polymerases I, II, and III subunit ABC3 Alternative name(s): DNA-directed RNA polymerase II subunit H DNA-directed RNA polymerases I, II, and I
<b>Expression Region</b>	2-150
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes one of the essential subunits of RNA polymerase II that is shared by the other two eukaryotic DNA-directed RNA polymerases, I and III.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.