



CAMSAP2 Antibody

Product Code	CSB-PA004478GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08AD1
Immunogen	Human CAMSAP1L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	calmodulin regulated spectrin-associated protein 1-like 1;CAMSAP1L1;KIAA1078;MGC150680;MGC150681 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CAMSAP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.