



CAP2 Antibody

Product Code	CSB-PA004487GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40123
Immunogen	Human CAP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	CAP, adenylate cyclase-associated protein, 2 (yeast);
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CAP2
Target Details	This gene was identified by its similarity to the gene for human adenylyl cyclase-associated protein. The function of This protein is unknown. However, the protein appears to be able to interact with adenylyl cyclase-associated protein and actin.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.