



CAPRIN2 Antibody

Product Code	CSB-PA004506GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IMN6
Immunogen	Human CAPRIN2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	caprin family member 2;CAPRIN2;C1QDC1;EEG-1;EEG1;FLJ11391;FLJ22569;KIAA1873;MGC10289 4;MGC134847;MGC134848;caprin-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CAPRIN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.