



CAPZA1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA004510GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P52907 |
| Immunogen | Human CAPZA1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | capping protein (actin filament) muscle Z-line, alpha 1;CAPZA1;CAPPA1;CAPZ;CAZ1 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | CAPZA1 |
| Target Details | CAPZA1 is a member of the F-actin capping protein alpha subunit family. This gene encodes the alpha subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |