



CASK Antibody

Product Code	CSB-PA004539GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14936
Immunogen	Human CASK
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	calcium/calmodulin-dependent serine protein kinase (MAGUK family);CASK;CAGH39;CMG;FGS4;FLJ22219;FLJ31914;LIN2;MICPCH;TNRC8 ;
Product Type	Purified Rabbit Anti Human,Mouse,Rat PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CASK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.