



CCDC115 Antibody

Product Code	CSB-PA004636GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NT0
Immunogen	Human CCDC115
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	coiled-coil domain containing 115;CCDC115;FLJ30131;MGC12981;ccp1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCDC115
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.