



# Recombinant *Xenopus tropicalis* Target of rapamycin complex 2 subunit MAPKAP1 (mapkap1), partial

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP013472XBF   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | A4IIM3  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | mapkap1   |
| <b>Protein Names</b>     | Recommended name: Target of rapamycin complex 2 subunit MAPKAP1 Short name= TORC2 subunit MAPKAP1 Alternative name(s): Mitogen-activated protein kinase 2-associated protein 1 Stress-activated map kinase-interacting protein 1 Short n  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Partial   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |