



CCT4 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA004859GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P50991 |
| Immunogen | Human CCT4 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | chaperonin containing TCP1, subunit 4 (delta);CCT4;CCT-DELTA;Cctd;MGC126164;MGC126165;SRB ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | CCT4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |