



CCT8 Antibody

Product Code	CSB-PA004873GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50990
Immunogen	Human CCT8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	chaperonin containing TCP1, subunit 8 (theta);CCT8;C21orf112;Cctq;D21S246;KIAA0002;PRED71 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCT8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.