



Recombinant *Geobacillus thermodenitrificans* Elongation factor Ts (tsf)

Product Code	CSB-EP025124GEH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4IMC4
Product Type	Recombinant Protein
Immunogen Species	<i>Geobacillus thermodenitrificans</i> (strain NG80-2)
Purity	>85% (SDS-PAGE)
Sequence	MAIT AQMVKE LREKTGAGMM DCKKALTETN GDMEKAIDWL REKGIKAAK KADRIAAEGM TYIATEGNAA VILEVNSETD FVAKNEAFQT LVKELAAHLL KQKPATLDEA LGQTMSSGST VQDYINEAVA KIGEKITLRR FAVVNKADDE TFGAYLHMGG RIGVLTLLAG SATEEVAKDV AMHIAALHPK YVSRDEVPQE EIAREREVLK QQALNEGKPE NIVEKMVEGR LKKFYEDVCL LEQAFVKNPD VTVRQYVESS GATVKQFIRY EVGEGLEKRQ DNFAEEVMSQ VRKQ
Source	E.coli
Target Names	tsf
Protein Names	Recommended name: Elongation factor Ts Short name= EF-Ts
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.