



CD164 Antibody

Product Code	CSB-PA004884GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04900
Immunogen	Human CD164
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	CD164 molecule, sialomucin;CD164;MGC-24;MUC-24;endolyn ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD164
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.