



CD84 Antibody

Product Code	CSB-PA004964GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIB8
Immunogen	Human CD84
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	CD84 molecule;CD84;DKFZp781E2378;LY9B;SLAMF5;hCD84;mCD84 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD84
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.